


LAMAGNETO PRO

LaMagneto Pro is the ideal PEMF device designed for those who have a family member with difficulty in using electromedical equipment, thanks to its range of pre-set programs that allow to start a therapy with just a simple click.



TECHNICAL FEATURES

- ✓ Two independent channels
- ✓ Colour display
- ✓ Magnetic field intensity up to 200 gauss per channel
- ✓ Frequency from 1 to 120 Hz
- ✓ 63 programs available
- ✓ Double user
- ✓ List of the last 10 programs used
- ✓ "Lock" function
- ✓ Counter of therapies performed and total hours of use

MEDICAL DEVICE CLASS IIa

LaMagneto Pro is a medical device CE1936

PROGRAMS

01 Osteoporosis	22 Cartilage lesion	43 Night 15 Hz
02 Arthrosis	23 Ligament lesion	44 Night 50 Hz
03 Arthritis	24 Bone oedema	45 Night 75 Hz
04 Cervical arthrosis	25 Tendinitis	46 Autoscan
05 Articular pain	26 Chondropathy	47 Treat. 1 Hz
06 Cervicalgias	27 Antinflammatory	48 Treat. 3 Hz
07 Chronic pain	28 Whiplash Syndrome	49 Treat. 5 Hz
08 Fractures	29 Healing	50 Treat. 8 Hz
09 Epicondylitis	30 Cutaneous ulcers	51 Treat. 10 Hz
10 Epitrocleitis	31 Discopathy	52 Treat. 15 Hz
11 Pseudoarthrosis	32 Myalgia	53 Treat. 16 Hz
12 Lumbalgy	33 Neuropathy	54 Treat. 20 Hz
13 Lumbar pain	34 Muscle strain	55 Treat. 30 Hz
14 Shoulder arthrosis	35 Muscle cramp	56 Treat. 40 Hz
15 Knee arthrosis	36 Rizarthrosis	57 Treat. 50 Hz
16 Scapulo Humeral Periarthritis	37 Impingement Syndrome	58 Treat. 60 Hz
17 Coxarthrosis	38 Carpal Tunnel	59 Treat. 70 Hz
18 Muscular atrophy	39 Titanium prosthesis	60 Treat. 75 Hz
19 Muscular treatment	40 Rotator cuff	61 Treat. 80 Hz
20 Osteonecrosis	41 Tarsal Tunnel	62 Treat. 90 Hz
21 Algodystrophy	42 Night 1 Hz	63 Treat. 100 Hz

LAMAGNETO PRO

KIT

- ✓ LaMagneto Pro
- ✓ A 3 solenoids belt applicator
- ✓ A professional applicator with a pair of solenoids
- ✓ Medical power supply
- ✓ User manual
- ✓ Magnet for verifying therapy operation
- ✓ Carrying bag
- ✓ Three non-woven fabric strips 15x150 cm

DIMENSIONS 180 x 110 x 50 mm

