

# T-ONE COACH

T-One Coach is a 4-channel electrotherapy device recommended for occasional or amateur athletes for trauma prevention and muscle strengthening.

## TECHNICAL FEATURES

- ✓ 4 connection cables (8 connectable electrodes)
- ✓ Constant current generator
- ✓ Compensated biphasic square wave
- ✓ 2+2 users programs
- ✓ Backlight graphic colour display
- ✓ 4,8 V rechargeable battery pack
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Frequency from 1 to 150 Hz
- ✓ Width impulse from 20 to 300  $\mu$ s
- ✓ Pre-adjusted time for each program



### MEDICAL DEVICE CLASS IIa

T-ONE COACH is a medical device CE0068

## KIT

- ✓ T-One Coach
- ✓ 4 connection cables (8 connectable electrodes)
- ✓ Self-adhesive square electrodes
- ✓ Self-adhesive rectangular electrodes
- ✓ Self-adhesive round electrodes
- ✓ Ni-MH rechargeable battery
- ✓ Battery charger
- ✓ User manual
- ✓ Electrodes position manual
- ✓ Carrying bag

## DIMENSIONS

150 x 90 x 35 mm

# T-ONE COACH

## PROGRAMS

| 16 REHA PROGRAMS      | 14 NEMS PROGRAMS               | 18 BEAUTY PROGRAMS                 | 2+2 USERS        |
|-----------------------|--------------------------------|------------------------------------|------------------|
| Tens rapid            | Warming up                     | Firming up                         | Tens rapid       |
| Tens endorphinic      | Resistance                     | Toning up                          | Tens endorphinic |
| Tens max values       | Resistance strength            | Definition                         | Firming up       |
| Inflammatory          | Basic strength                 | Modelling                          | Toning up        |
| Neck pain             | Fast strength                  | Microlifting                       | Definition       |
| Backache/Sciatic pain | Explosive strength             | Lipolysis abdomen                  | Modelling        |
| Sprains/Bruises       | Deep capillarization           | Lipolysis thighs                   | Microlifting     |
| Vascularization       | Muscle recovery                | Lipolysis glutei and hips          | Muscle relaxant  |
| Muscle relaxant       | Agonist/Antagonist             | Lipolysis arms                     | Massage          |
| Hand and wrist pain   | Sequential tonic contractions  | Tissue elasticity                  | Capillarization  |
| Plantar stimulation   | Sequential phasic contractions | Capillarization                    |                  |
| Epicondylitis         | Deep massage                   | Heaviness in legs                  |                  |
| Epitroclea            | Massaggio profondo             | Face Microlifting Handy Trode      |                  |
| Periarthritis         | EMS                            | Breast Microlifting Handy Trode    |                  |
| Tens Burst            |                                | Décolleté Microlifting Handy Trode |                  |
| Microcurrent          |                                | Face Wrinkles Handy Trode          |                  |
|                       |                                | Forehead Wrinkles Handy Trode      |                  |
|                       |                                | Collagen                           |                  |