

# Product datasheet

## G2 Duo

*The G2 Duo is the basis for welcoming your clients in a comfortable and safe manner. This treatment and examination couch can be moved fully downwards to chair height so that getting on or getting off has never been easier.*

*The long back section supports the sitting position, so that your client can be optimally involved in the treatment.*

### Standard features

- Two-section couch with adjustable back section
- Back section supported by gas spring
- Electric height adjustment with foot switch
- Cushions with 40 mm padding
- 15 cm clearance underneath frame for use with patient hoist (not compatible with 150 007)
- Non slip feet
- Cushions in a choice of 11 colours

### Ordering data

160 000 G2 Duo


### Optional

- 150 029 G2 Hand control
- 150 007 G2 Peripheral switch
- 150 001 G2 Hydraulic
- 150 062 G2 Wheel mechanism
- 150 003 G2 Cushion with face opening
- 150 004 G2 Cushions in non-standard colour

### Accessories

- 150 106 G2 Paper roll holder 50 cm
- 150 117 G2 Paper roll holder 60 cm
- 337 018 Paper rolls 50 cm (per 12 pieces)
- 332 915 Foot end protector for cushion
- xxx xxx Terry-cloth cover, stretch material (available in 9 colours 'Style')

### Specifications

Length x width:	± 205 x 67cm
Height adjustment:	± 47 – 100 cm
Weight:	± 85 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	230 VAC, 50-60 Hz, 1,0 A 120 VAC, 60 Hz, 1,7 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 MDD 93/42/EEC
	
Guarantee:	3 years (upholstery 1 year)




**G2 Duo**



Duo G2 shown with options:

- Peripheral switch
- Wheel mechanism
- Cushion with face opening

**Functionality/dimensions**



**Cushion colours**

**Style**



**Fashion**

